


Trimming and filtering



Filter data

© Can Stock Photo - csp34642023



Trim data

<https://www.dreamstime.com/stock-illustration-man-trimming-bush-garden-image64030511>


Overview of this talk

Removing poor quality data

Data decontamination

Synchronization of paired-end reads

Merging overlapping reads pairs



Why is it important to perform QC and filtering/trimming?

Data analysis also costs money and time

Filtering data

Length related

Quality score related

GC content related

Ambiguity code related

Sequence complexity related

Trimming data


Trim by length/position – fixed, e.g. 20 bp

Trim tails

Trim ends by quality scores


Illumina adapter sequences

When preparing a (TrueSeq) library, adaptors are ligated to the DNA of interest



We filter reads to remove poor quality data

Common to set some quality thresholds and remove all reads that does not comply
For example remove all reads with a lower average Q score than 20, or reads with more than 5 Ns




For example:
Remove all reads with a lower average Q score than 20

For example:
Remove all reads with more than 5 Ns

AfterQC - Automatic Filtering, Trimming, Error Removing and Quality Control for FASTQ data

Performs quality control and filtering/trimming of the sequence reads




AfterQC analyses the overlap of paired sequences for pair-end sequencing data

AfterQC will correct the low quality base according to its high quality mate




Removal of 3' adapter in the tail



Popular trimming tools

Many tools available – here are some :

- AfterQC
- Fastp
- Trimmomatic
- CutAdapt
- AlienTrimmer
- Sickle
- Trim Galore
- Sycthe
- Prinseq



It is important to remove sequence contaminations as early as possible

There can be many sources of contamination in the final sequence library

For example:
PhiX sequences from
the sequencing kit



For example:
Metagenomic samples may
contain sequences from the host

PhiX control

PhiX is used as a quality and calibration control for Illumina sequencing runs
10% of the genomes that are published in literature are contaminated with PhiX



<https://www.illumina.com>

Commentary | [Open Access](#) | [Published: 30 March 2015](#)

Large-scale contamination of microbial isolate genomes by Illumina PhiX control


[Supratim Mukherjee](#) , [Marcel Huntemann](#), [Natalia Ivanova](#), [Nikos C Kyrpides](#) & [Amrita Pati](#)

[Standards in Genomic Sciences](#) 10, Article number: 18 (2015) | [Cite this article](#)

FastQ Screen allows you to screen sequences against a set of sequence databases

Normally used to screen for host contaminations


Label sequences that match sequences in the database you provide



Popular decontamination tools


Removal of host contamination

- Fastq Screen
- DeconSeq
- CS-SCORE
- VecScreen
- Kraken




Synchronized FASTQ files are important for many tools

Some tools may remove one of the reads in a read-pair, hence the FASTQ files gets out of sync




Synchronized FASTQ files are important for many tools

Some tools may remove one of the reads in a read-pair, hence the FASTQ files gets out of sync




Synchronize FASTQ files with Repair

Repair pulls out the “singletons” and produce synchronized paired-end FASTQ files




DNA insert sizes and overlapping read pairs



Inner distance is the unsequenced middle piece


Insert size is stretch of sequence between the paired-end adapters

Fragment size includes the ligated Illumina adapters




When read pairs overlap, the inner distance is a negative number or a positive number describing the overlap

DNA insert sizes and overlapping read pairs



Inner distance = 80 bp
Insert size = $R_1 + R_2 + \text{inner distance} = 680$ bp
Fragment size = $R_1 + R_2 + \text{inner distance} + \text{adapters (x2)} = 800$ bp

Insert size in BBMerge is a bit confusing since the inner distance for overlapping reads is negative



Inner distance = BBMerge Insert size = -80
Insert size = R1 + R2 + BBMerge Insert size = 520 bp
Fragment size = R1 + R2 + inner distance + adapters (x2) = 640 bp

#InsertSize	Count (reads)
...	...
101	10248
102	10619
103	10397
104	10357
105	10218
106	10435
107	10288
108	9973
...	...

Generate longer reads by overlapping and merging read pairs before assembling a genome using BBmerge

Merging reads will improve the quality of the reads and generate longer reads

Longer reads allow the use of longer k-mers or fewer comparisons

Program	NA50 (bp)	Total Misassemblies	Indels/ 100 kbp	Genome Completeness (%)
Raw Data	60007	119	1.13	84.5
BBMerge	102577	127	0.84	84.88
BBMerge-REM	119328	117	0.81	85.18
BBMerge-RSEM	104441	115	0.84	84.88
COPE	89603	294	1.52	85.17
COPE-M3	98240	227	1.24	83.92
fastq-join	80672	183	1.17	84.74
FLASH	94846	282	1.41	85.20
leeHom	101992	290	1.1	84.91
PEAR	60937	660	1.46	84.28
Stitch	5623	20986	47.78	68.38
USEARCH	102156	131	0.88	84.77
XORRO	97403	158	1.08	84.85

Generate report and show to your boss 😊

MultiQC is a reporting tool that parses summary statistics from results and log files generated by other bioinformatics tools

Parses relevant information from log files to a HTML report file

The screenshot displays the MultiQC v1.3 web interface. On the left is a navigation sidebar with categories like 'General Stats', 'QUAST', 'Assembly Statistics', 'FastQC', and 'Adapter Content'. The main content area features the MultiQC logo, a description of the tool, and a 'General Statistics' table. The table lists various samples with their N50 (Kbp), Length (Mbp), and % Dups. A 'MultiQC Toolbox' is visible on the right side of the interface.

Sample Name	N50 (Kbp)	Length (Mbp)	% Dups
clean_megahit	3.4bp	29.6bp	
clean_metaspades	9.4bp	30.2bp	
sample_R1			0.0%
sample_R2			0.0%
sample_megahit	3.4bp	29.7bp	
sample_metaspades	4.0bp	29.7bp	
sample_trim_megahit	3.7bp	29.5bp	
sample_trim_metaspades	4.0bp	29.4bp	